

Investigating the Impact of Astrocyte-Melanoma Direct Junctions on Spatial GTP Distribution in Melanoma Cells Using GTP Biosensor

11/25/2024 BME milestone exam

Junqi Lu


Department of Biomedical Engineering

Mentor: Dr. Mikhail A. Nikiforov

Department of Biomedical Engineering

MBM exhibits high mortality & prevalence


- Melanoma represents **1%** skin cancer cases but **80%** skin cancer-related death due to high metastatic potential
- 50% metastatic melanoma patients develop brain metastasis (**MBM**)
- MBM patients' median survival = **12.8 month** with mortality rate = **80-85%**




MBM exhibits high mortality & prevalence

- Melanoma represents **1%** skin cancer cases but **80%** skin cancer-related death due to high metastatic potential
- 50% metastatic melanoma patients develop brain metastasis (**MBM**)



Understanding acquisition of invasive phenotype (melanoma → MBM) is critical for developing effective therapeutic strategies




Fluctuations in intracellular GTP levels regulate melanoma cell invasion & metastasis



Fluctuations in intracellular GTP levels regulate melanoma cell invasion & metastasis




GEVAL measures GTP distribution in melanoma cells




Anna Bianchi-Smiraglia,
Ph.D.

GEVAL measures GTP distribution in melanoma cells



GTP evaluators (GEVAL) structure
predicted by AlphaFold with confidence from **high** to **low**



GEVAL measures GTP distribution in melanoma cells



GEVAL measures GTP distribution in melanoma cells

GEVAL30



$K_{\text{eff}} = 32.3 \mu\text{M}$





GEVAL measures GTP distribution in melanoma cells

GEVAL30

$K_{\text{eff}} = 32.3 \mu\text{M}$




GEVALNull




**Aim: use computational tools & AI to streamline
GEVAL protocol & optimize analysis pipeline**


GEVAL*Iris* accelerates image processing




GEVALIris accelerates image processing



GEVALIris accelerates image processing



<https://youtu.be/aysFCivxEpw?t=1>




File folder (10)

- Heatmap_images
- Montage_images
- Processed_stacks
- Background_NaN
- Histogram_data
- Histogram_images
- ROI
- ROI_overlay
- Shape_data
- Background_data

Text Document (1)

mode.txt




Total 102,252 features

Flowcytometry verified GEVALNull has higher mean GFP intensities

Uninfected


GEVALNull

GEVAL30




GEVALNull population is brighter than GEVAL30

IMPDH is major GTP biosynthesis enzyme



MPA & AVN are IMPDH inhibitors






AVN treatment doesn't alter cell morphology

DMSO




LowAVN

HighAVN

Null




30



Color range normalized to [1E-4, 2.5]

Low AVN = 0.025 μ M
High AVN = 0.05 μ M₁₈


GEVALytics heatmap data analysis pipeline




AVN decreases GEVAL signals in GEVAL30

Pooled histograms

GEVALNull




GEVAL30



Low AVN = 0.025 μM
High AVN = 0.05 μM ₂₀

GEVALNull not responsive to AVN treatment while higher AVN significantly dropped GEVAL30 measurements


Pooled raincloud plot



Low AVN = 0.025 μM
High AVN = 0.05 μM ₂₁


Large variance may indicate suboptimal GEVAL performance

Pooled raincloud plot




Low AVN = 0.025 μM
High AVN = 0.05 μM ₂₂

Large variance may indicate suboptimal GEVAL performance




- What factor(s) contribute to large variance in GEVAL performance?
- How to reduce variance by screening out suboptimal GEVAL images more effectively—and even better, do it automatically?




Low AVN = 0.025 μM
High AVN = 0.05 μM ₂₃

GEVALytics AI workflow






Feature selection to reduce dimensionality





Feature selection to reduce dimensionality


102 features





Remove features with correlation coefficients ≥ 0.95



43 features





Random forest (RF) classifier training & best model



**Best RF classifier has 1000
decision trees**




Best RF classifier is highly reliable & less likely to overfit





Best RF classifier identifies laser settings as most important features contributing to large variance





Most important features to distinguish aren't highly correlated = no leftover data structure to model





GEVAL*Iris* + GEVALytics = automatic screening eyes with expertise

- RF classifier highly trustworthy & generalizes very well
 - Optimal GEVAL performance requires fine-tuning on laser settings
 - Automatically screens out suboptimal images for future GEVAL users, **anytime & anywhere**
- Potential future directions
 - Bigger image set to avoid curse of dimensionality
 - Correlation \neq causation so need controlled experiments to determine

Nikiforov Lab—Where Engineering Solutions Meet Scientific Curiosity

- **PI:** Mikhail A. Nikiforov, Ph.D.
- **Members**
 - Zhannan Han
 - Zhibo Yan, Ph.D.
 - Maja Beus, Ph.D.
 - Yihui Wang, Ph.D.
 - Junqi Lu
 - Loren Weidenhammer
 - Robert Noble

